

ScanPod

Non-Contact Temperature Check



For Low-Traffic Places Such As



Small and medium
-sized offices



Restaurants



Cafes



Pubs

Fever Screening Non-Contact Thermometer

- Designed for measuring human body temperature.
- AI Face Detection
Temperature detection range 30°C-45°C (86°F-113°F)
- Accuracy level within $\pm 0.3^{\circ}\text{C}$ ($\pm 0.54^{\circ}\text{F}$) without needing additional accessory. This level of accuracy can be used to screen elevated temperature which has to be confirmed by a secondary regular thermometer.
- Provides audio alarm when detecting persons with normal and elevated temperature.
- Detection setting: Mask Detection/Stranger Mode

Specifications:

Camera Type:	Infrared Thermometer
Spray System:	No
Display:	8 inch 800 x 1280 resolution
Body Dimension:	32.4 x 12.8 cm
Footprint:	35 x 30 cm (13.7" x 12")
Height:	155 cm (61")
Stand:	Floor and Desktop
Box Dimension: (1 unit)	112 x 32 x 15 cm (44" x 12.5" x 5.9")
	18 kg (39.6 lbs)



Support Safe Re-Open
of Businesses

WWW.NMSIMAGING.COM 301.526.7456
SALES@NMSIMAGING.COM

